Polinyl

Sutura chirurgica in poliammide monofilamento (non assorbibile) Sutura quirúrgica en monofilamento de poliamida (no absorbible)

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Descrizione

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sweden & martina

Sweden&Martina S.p.A.
Via Veneto 10 - 35020 Due Carrare, (PD) Italy
Tel. +39.049.91.24.300 - Fax. +39.049.91.24.290
www.sweden-martina.com

Polinyl

Sterile monofilament polyamide surgical suture (non absorbable) Chirurgischer Faden aus monofilem Polyamid (nicht resorbierbar).

ISTRUCTION FOR USE

Description

POLINYL suture is a sterile non-absorbable surgical suture device of synthetic origin, composed of polyamide 6 and polyamide 6/6, which derives from the polycondensation of hexamethylenediamide and adipic acid. **POLINYL** suture is hematin dyed black (C.I. 75290) and is available in a range of gauge sizes (UPS/EP) and lengths, attached to stainless steel needles of different sizes and shapes.

As for most of monofilament suture, **POLINYL** is indicated to suture infected and/or contaminated tissue; in addition, thanks to its lack of adherence to tissue, **POLINYL** is successfully used as a pull-out suture in post-operative course. The suture is not absorbed nor subject to degradation by the action of enzymes and it maintains its own tensile strength in "vivo". As with any foreign body, it may elicit a slight inflammatory reaction followed by encapsulation of the suture from fibrous connective tissue.

Indications

POLINYL suture is indicated for use in any procedure where the use of a nonabsorbable suture is required for suturing tissue or ligation.

Contraindications/Precautions/Recommendations

POLINYL sutures are to be used by qualified personnel, experienced in technical and surgical procedures where the use of nonabsorbable sutures is required. Care should be taken when handling the suture material to avoid any bending or flattening to the thread due to the application of surgical instruments such as forceps or needleholders. Greater knot security requires standard surgical technique of flat and square knots, to be applied with additional throws as may be indicated by the surgical circumstances and the surgeon's experience. Surgical sutures are disposable medical devices and must be disposed of as hospital waste material.

Adverse reaction

As any suture material, **POLINYL** is a foreign body and may cause inflammatory reaction.

Content

36 sterile envelopes.

Warnings

Keep in cool dry place, away from any direct heat.

Expiration date is shown on each single envelope and on the outer box.

Disposable.

Do not use if package is damaged.

Sterility

POLINYL sutures are sterilized by Gamma Ray.

Sterile unless package is opened or damaged.

Do not resterilize.

Reserved for the exclusive use by surgical medical personnel, in surgery or in hospital setting.

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