Resorb

Sutura chirurgica all'acido glicolico intrecciata e rivestita (assorbibile) Sutura quirúrgica a base de ácido glicólico, trenzada y revestida (absorbible)

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Descrizione

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ISTRUCTION FOR USE

Description

RESORB suture is a synthetic absorbable sterile surgical suture device, composed of braided multifilament of polyglycolic acid. The braid of natural white or dyed violet (D&C Violet N°2 – C.I. 60725) multifilament is coated with a mixture of polycaprolactone and calcium stearate or poly(glycolide-co-L-lactide) and calcium stearate. **RESORB** sutures are available in a range of gauges sizes (USP/EP) and lengths, attached to stainless steel needles of different size and shapes. The suture **RESORB** can accidentally cause an initial mild inflammatory reaction in the tissues resulting in the formation of fibrous connective tissue. The progressive loss of the tensile strength of the suture occurs as a result of a hydrolysis process that degrades the polyglycolic acid, allowing it to be absorbed and subsequently metabolized by the body. Studies performed on animals show that the absorption of the **RESORB** suture is practically complete around the 90th day after implantation (medium term absorption means absorption after 60-90 days)

Post implantation	Tensile strength residual
1 week	98%
2 weeks	65%
3 weeks	35%
4 weeks	5%

Further details are described in the general catalog. **RESORB** suture meets the requirements of the Pharmacopoeia European Union for Braided Sterile Absorbable Synthetic Sutures.

Indications

RESORB suture is indicated in the approximation of soft tissues and in all interventions that at the discretion of the user require an absorbable suture.

Contraindications

As all the absorbable sutures, **RESORB** should not be used when long-term wound support is required.

Precautions

RESORB sutures are to be used in hospitals operating theatres and ambulatories or in any hospital environment by qualified personnel, experienced in technical and surgical procedures where the use of an absorbable suture is required. It is therefore recommended that the surgeon will consider the possible use of non absorbable sutures for the closure of tissues which may undergo distension or stretching or requiring additional support. This suture kind may be inappropriate in elderly, malnourished or weakened patients or in patients suffering from conditions which may delay wound healing process.

Recommendations

Care should be taken when handling the suture material to avoid any bending or flattening to the thread due to the application of surgical instruments such as forceps or needleholders. Greater knot security requires standard surgical technique of flat and square knots, to be applied with additional throws as may be indicated by the surgical circumstances and the surgeon's experience. Surgical sutures are disposable medical device and must be disposed of as hospital waste material.

Adverse reaction

RESORB suture may cause local inflammatory reaction associated to erythema formation.

Content

36 sterile envelopes.

Warnings

Keep in cool dry place, away from any direct heat.

Expiration date is shown on each single envelope and on the outer box.

Disposable.

Do not use if package is damaged.

Sterility

RESORB sutures are sterilized by ethylene oxide.

Sterile unless package is opened or damaged.

Do not resterilize.

Reserved for the exclusive use by surgical medical personnel, in surgery or in hospital setting.

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Sweden & Martina S.p.A.