NUMERI UNO

ESSE&EMME NEWS MAGAZINE

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Fatigue resistance of nickel-titanium rotary instruments in continuous or reciprocating rotation

Introduction

To investigate cyclic fatigue resistance of four different endodontic rotary instruments produced by M-wire nickel - titanium (Ni-Ti), R-Phase Ni-Ti or traditional Ni-Ti tested in continuous or reciprocating rotation.

Materials and Methods

A total of 180 new Reciproc R25; WaveOne Primary; Mtwo size 25.06; Twisted File 25.06 were tested for cyclic fatigue resistance. Forty-five files of the same brand were randomly assigned to three different groups of 15 each. Instruments were tested using a specific reciprocating torque-controlled motor*. Files of group 1 were tested in continuous rotation; group 2 in reciprocating motion "Reciproc ALL"; group 3 in reciprocating motion "WaveOne ALL". Continuous rotation was performed at 300 rpm as in clinical practice, while reciprocating motions "Reciproc ALL" and "WaveOne ALL" are preset, unchangeable and patented by the manufacturer. Cyclic fatigue tests were performed by a mechanical device already used in previous studies. The apparatus enabled the instrument to rotate freely within a stainless steel artificial canal (60 degree angle of curvature and 5 mm radius of curvature) manufactured by reproducing instrument's size and taper. For each files, numbers of cycles to failure (NCF) was recorded. A two-way analysis of variance and Bonferroni multiple comparison post-hoc tests at 0.05 were used to evaluate the data.

Results

The type of rotation with pre-programmed reciprocating movement "Reciproc ALL" or "WaveOne ALL" - group 2 or 3 - showed mean of NCF values higher than that one of continuous rotation group for the same brand (Two way ANOVA $P < 0.0001$; interaction 0,08). Significant difference was not found between the NCF of the two reciprocating rotation tested (group 2 and 3) for all instruments ($P > 0.05$). The NCF mean of all Mtwo and Twisted File groups is greater than that of Reciproc ($P < 0.05$ for each comparison) or WaveOne ($P < 0.0001$ for each comparison). The NCF was not significant different between Mtwo or Twisted File groups. Cyclic fatigue resistance of Reciproc groups was statistically significant higher than the same groups of WaveOne ($P < 0.001$) except for group 2 - files used with pre-programmed reciprocating movement "Reciproc ALL" - ($P > 0.05$).

Conclusions

Reciprocating movements enhance the NCF of all instruments tested. The pre-programmed reciprocating movement "Reciproc ALL" increases cyclic fatigue resistance of files significantly. In these experimental conditions, Twisted File 25.06 and Mtwo 25.06 had a higher resistance to cyclic fatigue with any type of rotation than Reciproc R25 and WaveOne Primary.

* VDW Silver Reciproc, VDW GmbH – Dentsply International Inc., Munich, Germany

Key-words: Cyclic fatigue, M-wire, R-Phase, Nickel - titanium; reciprocating rotation.

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Influence of operational speed on cyclic fatigue resistance of Mtwo nickel - titanium instruments

Introduction

To evaluate the influence of different operational speeds on cyclic fatigue resistance of Mtwo Nickel - titanium (Ni-Ti) instruments.

Methods and Materials

A total of 180 new Mtwo Ni-Ti rotary files different in size and taper were tested. Forty-five instruments for each type (10.04, 15.05, 20.06 and 25.06) were randomly assigned to three different groups of 15 each. Instruments were used in the sequence suggested by the manufacturer to shaping 9 standardized simulated canals at 350 rpm (group 1); at 250 rpm (group 2); at 150 rpm (group 3). The time (seconds) needed to shaping 9 simulated canals for every instrument was recorded.

After this shaping phase, cyclic fatigue resistance of all instruments were tested. Cyclic fatigue tests were performed by a mechanical device already used in previous studies. The apparatus enabled the instrument to rotate freely (at the same speed used in the shaping phase) within a stainless steel artificial canal (60 degree angle of curvature and 5 mm

radius of curvature) manufactured by reproducing instrument's size and taper. For each files, time to fracture (TtF) from the start of the test until the moment of file breakage was recorded. A two-way analysis of variance and Bonferroni multiple comparison post-hoc tests at 0.05 were used to evaluate the data.

Results

Significant difference was found in the TtF according to size and taper of the files. Instruments 10.04 showed the highest cyclic fatigue resistance ($P < 0.0001$ vs 25.06 files). The time needs to shaping 9 simulated canals was significantly lower at 350 rpm than at 150 rpm. The cyclic fatigue resistance wasn't significantly different ($P > 0.05$) between the three groups of the same file's type considering the operational speed as independent variable.

Conclusions

Cyclic fatigue resistance was influenced by the type of file (size and taper), but it was not related to the speed at which the NiTi files were operated.

Key-words: *Cyclic fatigue; nickel-titanium, operational speed; rotary instruments; simulated canals.*

Influence of fiber post on fracture resistance of molar restored with composite onlay

Introduction

Previous studies have emphasized that endodontic treatment is a major etiological factor for tooth fractures.

The purpose of this in vitro study was to evaluate the effect of post placement on the fracture resistance of endodontically treated permanent teeth. The null hypothesis was that there would be significant differences among experimental groups restored with and without fiber posts in indirect overlay composite restoration.

Materials and Methods

Fifty noncarious mandibular molars with mature apices, extracted for periodontal reasons, were selected. Endodontic treatment was carried out in all samples excluding the control group (intact teeth, group 5).

After 24 h, the samples were randomly divided into five groups (n = 10 each). In groups 1 and 2, cavities with two walls remaining were prepared. In groups 3 and 4, cavities with only the buccal wall remaining were prepared. Thus, the residual thickness of the remaining walls at the height of the contour was 2 ± 0.2 mm. In groups 1 and 3, samples were treated with All-Bond 3 (Bisco, Schaumburg, IL, USA) and a build-up was performed with A2-shaded nanohybrid resin composite (Venus Diamond; Heraeus Kulzer) using an oblique layering technique.

In groups 2 and 4 a post space was prepared in the distal canal of each sample to a depth of 5 mm using a dedicated drill (Dentsply Maillefer). Adhesive procedures were performed with All-Bond 3 (Bisco) as described for groups 1 and 3. A dual-curing cement (Core X flow; Dentsply Caulk, Milford, DE, USA) was placed into the canal. Number 2 Radix fiber posts (Dentsply Maillefer) were cemented. A composite build-up restoration was then performed as described in groups 1 and 3.

Standardized adhesive overlay preparation was performed for all samples, reducing the remaining cusps and composite build-up to 2 mm in all directions.

The overlays were cemented using a hydrophobic light-curing bonding resin (OptiBond FL; Kerr) and then a light-curing resin composite (Venus Diamond; Heraeus Kulzer) was inserted into all cavities.

All specimens were subjected to 3000 thermal cycles between 5°C and 55°C for 60s each.

All groups were exposed to cyclic loading (Mini Bionix II; MTS Systems, Eden Prairie, MN, USA) with an inclination angle of 45° to the long axis of the tooth at a frequency of 8 Hz, start-

ing with a load of 20 N for 5000 cycles, followed by a load of 50 N at a maximum of 20,000 cycles. A 2-mm-diameter metallic ball was used. The site of loading was the central fossa. The specimens were then submitted to the static fracture resistance test using a universal testing machine (Instron, Canton, MA, USA) with a 2-mm-diameter steel sphere crosshead welded to a tapered shaft and applied to the specimens at a constant speed of 2 mm/min and an angle of 45° to the long axis of the tooth. The forces necessary to fracture each tooth were measured in Newtons (N). Fractured specimens were observed under optical microscopy (Wild, Heerbrugg, Gaiss, Switzerland) at 40X magnification to establish failure modes, which were classified as restorable failures, including adhesive failures above the CEJ; or non-restorable failures, including vertical root fractures below the CEJ.

Data are expressed as means \pm standard deviations and frequencies (%). The Kolmogorov-Smirnov statistical test for normality revealed a normal data distribution. The Kruskal-Wallis test was used for comparison among groups, and the χ^2 test was used to compare the failure modes of the specimens. Differences were considered statistically significant when $p < 0.05$. All statistical analyses were performed using the SPSS software (ver. 19.0 for Windows; SPSS Inc., Chicago, IL, USA).

Results

The mean values of fracture resistance in all groups are listed in Table 1. The Kruskal-Wallis test showed no significant difference among groups ($H = 7.814$, $p = 0.065$). The evaluation of failure mode among all groups is reported in Table 1. The χ^2 test showed no significant difference in the fracture patterns among groups ($\chi^2 = 1.6368$, $p = 0.200$). The highest proportion of restorable fractures (70%) was observed in groups 2 and 4, and the highest proportion of non-restorable fractures (60%) was observed in group 1.

Discussion

In the current study, we examined the fracture resistance of endodontically treated mandibular molars restored with composite overlays supported by build-up with or without fiber post insertion. The data led to the conclusion that the initial null hypothesis must be rejected because the fracture values of different restorations tested did not differ significantly.

The results of this in vitro study show that the insertion of a glass fiber-reinforced post during the build up of mandibu-

lar molars with one or two residual cavity walls and overlay preparation margins coronal to the CEJ appeared not to significantly increase the fracture resistance of the restoration.

Conclusions

Within the limits of this in vitro study, we can assert that fiber post placement does not fundamentally increase fracture

resistance of endodontically treated mandibular molars with indirect composite overlay restorations. However, fiber post insertion could be more helpful in distributing cervical or radicular stresses due to extra axial loads. Long-term clinical studies are required to confirm the present findings.

References

1. Lagouvardos P., Sourai P., Douvitsas G.; *Coronal fractures in posterior teeth*; Operative Dentistry 1989; 14(1): 28-32
2. Hürmüzlü F., Kiremitçi A., Serper A., Altundaşar E., Siso S.H.; *Fracture resistance of endodontically treated premolars restored with ormocer and packable composite*; Journal of Endodontics 2003; 29(12): 838-840
3. Hansen E.K., Asmussen E.; *In vivo fractures of endodontically treated posterior teeth restored with enamel-bobbed resin*; Endodontics and Dental Traumatology 1990; 6: 218-225
4. Assif D., Gorfil C.; *Biomechanical considerations in restoring endodontically treated teeth*; Journal of Prosthetic Dentistry 1994; 71(6): 565-567.
5. Kishen A.; *Mechanisms and risk factors for fracture predilection in endodontically treated teeth*; Endodontic Topics 2006; 13: 57-83
6. Tang W., Wu Y., Smales R.J.; *Identifying and reducing risks for potential fractures in endodontically treated teeth*; Journal of Endodontics 2010; 36(4): 609-617
7. Caplan D.J., Kolker J., Rivera E.M., Walton R.E.; *Relationship between number of proximal contacts and survival of root canal treated teeth*; International Endodontic Journal 2002; 35(2): 193-199.
8. Nagasiri R., Chitmongkolsuk S.; *Long-term survival of endodontically treated molars without crown coverage: a retrospective cohort study*; Journal of Prosthetic Dentistry 2005; 93(2): 164-170
9. Boschian Pest L, Guidotti S., Pietrabissa R., Gagliani M.; *Stress distribution in a post-restored tooth using the three-dimensional finite element method*; Journal of Oral Rehabilitation 2006; 33(9): 690-697
10. Ferrari M., Cagidiaco M.C., Grandini S., De Sanctis M., Goracci G.; *Post placement affects survival of endodontically Treated Premolars*; Journal of Dental Research 2007; 86(8): 729-734
11. Scotti N., Scansetti M., Rota R., Pera F., Pasqualini D., Berutti E.; *The effect of the post length and cusp coverage on the cycling and static load of endodontically treated maxillary premolars*; Clinical Oral Investigations 2011; 15(6): 923-929
12. Santos A.F., Meira J.B., Tanaka C.B., Xavier T.A., Ballester R.Y., Lima R.G., Pfeifer C.S., Versluis A.; *Can fiber posts increase root stresses and reduce fracture?*; Journal of Dental Research 2010; 89(6): 587-591
13. Mangold J.T., Kern M.; *Influence of glass-fiber posts on the fracture resistance and failure pattern of endodontically treated premolars with varying substance loss: an in vitro study*; Journal of Prosthetic Dentistry 2011; 105(6): 387-393
14. Newman M.P., Yaman P., Dennison J., Raftar M., Billy E.; *Fracture resistance of endodontically treated teeth restored with composite posts*; Journal of Prosthetic Dentistry 2003; 89(4): 360-367
15. Salameh Z., Sorrentino R., Papacchini F., Ounsi H.F., Tashkandi E., Goracci C., Ferrari M.; *Fracture resistance and failure patterns of endodontically treated mandibular molars restored using resin composite with or without translucent glass fiber posts*; Journal of Endodontics 2006; 32(8):752-755

Evaluation of antibacterial effects of Er, Cr: YSGG laser at different irradiation time on *Enterococcus faecalis* in infected root canals

Introduction

The success of endodontic treatment depends on the effective elimination of microorganisms from the root canal, and lasers demonstrated to provide effective disinfection when used as adjunctive device to the conventional treatment. The aim of this in vitro study was to determine the effectiveness of the Er, Cr: YSGG laser by measuring its bactericidal effect inside root canal experimentally colonized with *Enterococcus faecalis*. We also tested its antimicrobial effect at different irradiation time (30 and 60 seconds).

Materials and Methods

A total of 39 single-rooted extracted human teeth were endodontically prepared with rotary instrumentation. All were sterilized and inoculated with a suspension of *Enterococcus faecalis* (105 bacteria/ml). The teeth were randomized into two treatment (Group 1, Group 2) and one control groups. In all groups, teeth were chemical irrigated with 5.25% sodium hypochlorite and 17% ethylenediaminetetracetic acid (EDTA). Group 1 and 2 were also irradiated at 30 and 60 seconds respectively, with an Er, Cr: YSGG laser at 75 mJ and 10 Hz with a pulse duration of 140 μ s. The tip (diameter = 200 μ m) was placed into the root canal access and was kept stationary and

not advanced into the canal. The teeth were finally rinsed with saline water for 30 seconds. The teeth of control group were rinsed with saline water after 30 minutes, as recommended by the standard protocol. Paper points used to sample bacteria from the root canals were transferred to tubes containing brain heart infusion (BHI) broth. Tubes were incubated and the appearance of broth turbidity was indicative of bacteria remaining in the root canal. Microbiological procedures were performed in blindness.

Results

For all groups, a bactericidal effect was observed. The use of laser at 30 and 60 seconds eliminated a percentage of 92.3 and 100% of *Enterococcus faecalis*, respectively. In the control group a reduction of 92.3% was observed. No statistical differences were observed between the three groups ($P=0.543$, Fischer exact test).

Discussion

The results indicated a bactericidal effect of Er, Cr: YSGG laser irradiation at the settings used in this study. Moreover, the highest bactericidal effect of this laser was observed at 60 seconds irradiation time.

References

1. Nair P.N. (2004); *Pathogenesis of apical periodontitis and the causes of endodontic failures*; Crit Rev Oral Biol Med 15:348-381
2. Sundqvist G. (1998); *Microbiologic analysis of teeth with failed endodontic treatment and the outcome of conservative re-treatment*; Oral Surg Oral Med Oral Pathol Oral Radiol Endod 85:86-93
3. Distel J.W., Hatton J.F., Gillespie M.J. (2002); *Biofilm formation in medicated root canals*; J Endod 28:689-693
4. Klinke T., Klimm W., Gutknecht N. (1997); *Antibacterial effects of Nd:YAG laser irradiation within root canal dentine*; J Clin Laser Med Surg 15:29-31
5. Odor T.M., Chandler N.P., Watson T.F., Ford T.R., McDonald F. (1999); *Laser light transmission in teeth: a study of the patterns in different species*; Int Endod J 32:296-302
6. White J.M., Goodis H.E., Cohen J.N. (1991); *Bacterial reduction of contaminated dentine by Nd:YAG laser*; J Dent Res 70:412 [abstract 1170]
7. Gutknecht N., Moritz A., Conrads G., Sievert T., Lampert F. (1996); *Bactericidal effect of the Nd:YAG laser in vitro root canals*; J Clin Laser Med Surg 14:77-80
8. Moritz A., Doertbudak O., Gutknecht N., Goharkhay K., Schoop U., Sperr W. (1997); *Nd:YAG Laser irradiation of infected root canals in combination with microbiologic examinations-an in vivo study*; J Am Dent Assoc 128:1525-1530
9. Schoop U., Kluger W., Moritz A., Nedjelik N., Georgopoulos A., Sperr W. (2004); *Bactericidal effect of different laser systems in the deep layers of dentin*; Lasers Surg Med 35:111-116
10. Moritz A., Gutknecht N., Goharkhay K., Schoop U., Wernisch J., Sperr W. (1997); *In vitro irradiation of infected root canals with a diode laser: results of microbiologic, infrared spectrometric, and stain penetration examinations*; Quintessence Int 28/3:205-209
11. Schoop U., Kluger W., Dervisbegovic S., Goharkhay K., Wernisch J., Georgopoulos A., Sperr W., Moritz A. (2006); *Innovative wavelengths in endodontic treatment*; Lasers Surg Med 38:624-630
12. Schoop U., Moritz A., Kluger W., Patruta S., Goharkhay K., Sperr W., Wernisch J., Gattringer R., Mrass P., Georgopoulos A. (2002); *The Er:YAG laser in endodontics: results of an in vitro study*; Lasers Surg Med 30:360-364
13. Eversole L.R., Rizoiu I.M. (1995); *Preliminary investigations on the utility of an erbium, chromium YSGG laser*; J Calif Dent Assoc 23:41-47
14. Yamazaki R., Goya C., Yu D.G., Kimura Y., Matsumoto K. (2001); *Effects of erbium, chromium:YSGG laser irradiation on root canal walls: a scanning electron microscopic and thermographic study*; J Endod 27:9-12
15. Kimura Y., Yu D.G., Kinoshita J., Hossain M., Yokoyama K., Murakam Y., Nomura K., Takamura R., Matsumoto K. (2001); *Effects of erbium, chromium: YSGG laser irradiation on root surface: morphological and atomic analytical studies*; J Clin Laser Med Surg 19:69-72

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Void formation between carrier-based obturator core and root canal sealer due to stripped filling material

Introduction

The final step in endodontic therapy is root canal system obturation with the objective to achieve apical and coronal seal, sealing of accessory and lateral canals; all pathways from the root canal system to the surrounding periodontal ligament should be filled. Gutta percha combined with a root canal sealer cement has been the root canal system obturation of choice for many years. Carrier-based obturation materials provide a vehicle for delivery of the root canal filling material in one step. These materials are either delivered cold or thermosoftened. Conventional carrier-based systems rely on gutta percha as the filling material. The main problem of these obturation systems is stripping of the filling material from the carrier during the insertion of the obturator mainly due to constrictions in the canal. The aim of our study is to investigate whether the resin-based core of innovative carrier based obturators which use resin-based Resilon as filling material (*Real Seal 1, Sybron Dental Specialties, Orange, CA, USA*) can bond directly to a resin-based sealer.

Materials and Methods

34 patients who presented a tooth to be extracted for orthodontic or periodontal reasons were treated endodontically. All endodontic treatments were performed by the same operator and were completed in a single visit. In each case the operative field was isolated with rubber dam and an endodontic access cavity was prepared using a diamond truncated cone bur and an endo zeta tungsten non-end cutting bur mounted on a high-speed handpiece under abundant irrigation. Stainless steel hand files were used to survey the root canal system. Working length was determined by the joint use of both an electronic apex locator (*PROPEX II-DENTSPLY MAILLEFER, Ballaigues, CH*) and an intraoral intra-operative radiograph. Cleaning and shaping of the root canal system of each molar was performed using the simultaneous technique with Mtwo NI-TI rotary files (*Sweden & Martina SPA, Padova, IT*) and Mtwo Apical NI-TI rotary files (*Sweden & Martina SPA, Padova, IT*) for the preparation of the apical area. After each instrument the root canal system was irrigated alternating 5% sodium hypochlorite (NaOCl) and 17% EDTA. In cases of retreatment, Mtwo R 25/05 file (*Sweden & Martina SPA, Padova, IT*) was used in order to remove filling materials. After chemo-mechanical preparation the obturation of the root canal system of 17 teeth was carried out by

using a Real Seal 1 (*Sybron Dental Specialties, Orange, CA, USA*) obturator which prior to obturation was stripped off its Resilon cover layer. The rest 17 teeth were obturated using a "naked" Thermafil obturator (*DENTSPLY MAILLEFER, Ballaigues, CH*). A resin-based dual cure self-etch root canal sealer (*RealSeal SE, Sybron Dental Specialties, Orange, CA, USA*) was used together with all "naked" obturators of both groups during the obturation of the root canal systems. One week after endodontic treatment the teeth were extracted. Extracted teeth were prepared adequately examined under microscope ("Sapienza" University of Rome, Italy). The Statistical Package for the Social Sciences (*SPSS Inc., Chicago, Ill.*) was used to analyze the data.

Results

Traditional plastic Thermafill obturator core (*DENTSPLY MAILLEFER, Ballaigues, CH*) showed gaps and fissures between its surface and the root canal sealer (*RealSeal SE, Sybron Dental Specialties, Orange, CA, USA*).

On the other hand the surface of the resin-based Real Seal 1 core showed a seamless interface between its surface and the resin-based sealer (*RealSeal SE, Sybron Dental Specialties, Orange, CA, USA*) which bonds to the tubules and the canal dentin.

Discussions

The introduction of bonded obturator materials (methacrylate resins) has enabled the clinician to obtain a bonded seal to the root canal dentin in areas reached by the adhesive material. If the obturating material is gutta percha such bonding cannot occur. Gutta percha as a root canal filling material offers numerous advantages such as biocompatibility, compatibility with sealer materials, ability to be heated and softened for canal placement and technique flexibility, is inexpensive, has a long history of use and can be easily removed in cases of retreatment. Nonetheless doesn't provide a seal irrespective of the placement technique and must always used with a sealer. Sealers used with gutta percha (calcium hydroxide, calcium phosphate, zinc oxide eugenol, epoxy resin, methacrylate resin) fill voids between gutta percha points when multiple points are used and are responsible for obtaining apical and coronal seal.

Furthermore gutta percha if it has been heated prior to placement it can shrink as it cools. Moreover regarding carrier-

based obturation technique if gutta percha is stripped off a carrier based plastic core obturator during the obturation of a root canal system a gap is created that can compromise apical/coronal seal and furthermore the success of the whole endodontic treatment as bacteria can travel down this void to the apex. It is worth underlining that a new carrier-based obturator called Guttacore (*DENTSPLY Tulsa Dental Specialties, Tulsa, USA*) has been recently launched in the USA market; the difference of these obturators is that the core of the carrier is made out of gutta percha, as the filling material too, so it facilitates the clinician in the preparation of the post space and in cases of retreatment but it doesn't resolve the bonding question of the filling material to the root canal sealer and of course to root canal dentin.

Conclusions

Methacrylate-based obturator has resulted in the availability of a carrier-based obturation material that utilizes adhesive technology for obturation of the root canal system that can guarantee a monoblock bonded obturation and therefore apical and coronal sealing even if the filling material is stripped off during the obturation "correcting" in this way the main sideeffect of the obturation technique itself. Only a true chemical bond between root dentin, root canal sealer and root canal filling material could create a real seal preventing leakage and reinfection. Authors retain that future studies of great interest that could be conducted; concern fluid filtration leakage of the interface resin-based Real Seal 1 core (*Sybron Dental Specialties, Orange, CA, USA*) and resin-based sealer (*RealSeal SE, Sybron Dental Specialties, Orange, CA, USA*).

References

1. Mounce R., Glassman G.; *Bonded Endodontic Obturation: Another Quantum Leap Forward for Endodontics*; Oral Health; 2004 Jul; 13-22.
2. Epley S.R., Fleischman J., Hartwell G., Cicalese C.; *Completeness of root canal obturations: Epiphany techniques versus gutta-percha techniques*; J Endod. 2006 Jun;32(6): 541-4.
3. Verissimo D.M., do Vale M.S., Monteiro A.J.; *Comparison of apical leakage between canals filled with gutta percha/ AH-Plus and the Resilon/ Epiphany System, when submitted to two filling techniques*; J Endod 2007 Mar; 33(3): 291-4.
4. Ungor M., Onay E.O., Orucoglu H.; *Push-out bond strengths: the Epiphany-Resilon endodontic obturation system compared with different pairings of Epiphany, Resilon, AH Plus and gutta percha*; Int Endod J. 2006 Aug; 39 (8) : 647-7.
5. Vail et al.; *An in vitro Comparison of Microleakage between Resilon and Gutta Percha with a Fluid Filtration Model*; J Endod 2007 Dec.,Vol 33(12).

Importance of immediate post-endodontic coronal seal

Introduction

Analyse and evaluate the various strategies regarding post-endodontic restorative treatment.

Materials and Methods

280 adult patients who presented a mandibular molar to be treated endodontically were divided in 4 groups. All endodontic treatments were performed by the same operator and were completed in a single visit. In each case the operative field was isolated with rubber dam and an endodontic access cavity was prepared using a diamond truncated cone bur and an endo zeta tungsten non-end cutting bur mounted on a high-speed handpiece under abundant irrigation. Stainless steel hand files were used to survey the root canal system. Working length was determined by the joint use of both an electronic apex locator (*PROPEX II-DENTSPLY MAILLEFER, Ballaigues, CH*) and an intraoral intraoperative radiograph. Cleaning and shaping of the root canal system of each molar was performed using the simultaneous technique with Mtwo NI-TI rotary files (*Sweden & Martina SPA, Padova, IT*) and Mtwo Apical NI-TI rotary files (*Sweden & Martina SPA, Padova, IT*) for the preparation of the apical area. After each instrument the root canal system was irrigated alternating 5% sodium hypochlorite (NaOCl) and 2.5% EDTA. In cases of retreatment, Mtwo R 25/05 file (*Sweden & Martina SPA, Padova, IT*) was used in order to remove filling materials. After chemo-mechanical preparation the obturation of the root canal system of all teeth was carried out by using an obturation system in which the gutta-percha is preapplied to a carrier and is heated in a special heater (*Domino, Sweden & Martina SPA, Padova, IT*). After radiographic evaluation of the root canal system obturation a different post-endodontic restorative treatment was performed to each group. Group I: Immediate application of prefabricated glass fiber post (*Mtwo KOR, Sweden & Martina SPA, Padova, IT + VI-RAGE DUAL PLUS, Sweden & Martina SPA, Padova, IT*) and coronal composite adhesive restoration (*ACE ALL BOND SE, BISCO DENTAL PRODUCTS CANADA INC, BC, Canada + BIS CORE, BISCO DENTAL PRODUCTS CANADA INC,*

BC, Canada + ADONIS, Sweden & Martina SPA, Padova, IT). Group II: Immediate coronal composite adhesive restoration (*ACE ALL BOND SE, BISCO DENTAL PRODUCTS CANADA INC, BC, Canada + ADONIS, Sweden & Martina SPA, Padova, IT*). Group III: Application of an adhesive system (*ACE ALL BOND SE, BISCO DENTAL PRODUCTS CANADA INC, BC, Canada*) and a temporary filling. Coronal composite adhesive restoration (*ACE ALL BOND SE, BISCO DENTAL PRODUCTS CANADA INC, BC, Canada + ADONIS, Sweden & Martina SPA, Padova, IT*) performed in a second visit a week after the root canal system obturation. Group IV: Application of a temporary filling (*Cavit, 3M ESPE, USA*). Coronal composite adhesive restoration (*ACE ALL BOND SE, BISCO DENTAL PRODUCTS CANADA INC, BC, Canada + ADONIS, Sweden & Martina SPA, Padova, IT*) performed in a second visit a week after the root canal system obturation. The Statistical Package for the Social Sciences (*SPSS Inc., Chicago, Ill.*) was used to analyze the data.

Results

Statistically significant differences regarding the success rate of teeth occurred between the first three groups and the fourth one.

Discussions

Temporary filling alone cannot guarantee a coronal sealing after an endodontic treatment. Authors retain in accordance with international literature that it would be ideal after an endodontic treatment to proceed with an immediate composite adhesive restoration with or without an application of a prefabricated glass fiber post. Single visit endodontic treatments can be followed by an immediate coronal adhesive coronal restoration thanks to the innovative endodontic and adhesive restoration materials that minimize operative time.

Conclusions

A non immediate post endodontic access cavity sealing can compromise the success of a valid (homogeneous filling material throughout the entire length of each root canal reaching radiographic apex) root canal obturation.

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Bibliografia

1. Amaral M., Rippe M.P., Bergoli C.D., Monaco C., Valandro L.F.; *Multi-step adhesive cementation versus one-step adhesive cementation: push-out bond strength between fiber post and root dentin before and after mechanical cycling*; Faculty of Odontology, Sao Paulo State University (UNESP), Sao Jose dos Campos, Brazil; Gen Dent. 2011 Sep-Oct.
2. Stawarczyk B., Basler T., Ender A., Roos M., Ozcan M., Hämmerle C.; *Effect of surface conditioning with airborne-particle abrasion on the tensile strength of polymeric CAD/CAM crowns luted with self-adhesive and conventional resin cements*; Clinic of Fixed and Removable Prosthodontics and Dental Material Science, Center of Dental Medicine, University of Zurich, Zurich, Switzerland; J Prosthet Dent. 2012 Feb
3. Cortellini D, Canale A.; *Bonding lithium disilicate ceramic to feather-edge tooth preparations: a minimally invasive treatment concept*; J Adhes Dent. 2012 Feb
4. Beier U.S., Kapferer I., Dumfahrt H.; *Clinical long-term evaluation and failure characteristics of 1,335 all-ceramic restorations*; Clinical Department of Restorative and Prosthetic Dentistry, Innsbruck Medical University, Austria. Int J Prosthodont. 2012 Jan-Feb
5. Taschner M., Krämer N., Lohbauer U., Pelka M., Breschi L., Petschelt A., Frankenberger R.; *Leucite-reinforced glass ceramic inlays luted with self-adhesive resin cement: A 2-year in vivo study*; Dental Clinic 1 - Operative Dentistry and Periodontology, University Medical Center Erlangen, University of Erlangen-Nuremberg Dent Mater. 2012 May
6. Gresnigt M., Ozcan M.; *Esthetic rehabilitation of anterior teeth with porcelain laminates and sectional veneers*; University Medical Center Groningen, Center For Dentistry and Dental Hygiene, Department of Fixed and Removable Prosthodontics, Groningen, The Netherlands. J Can Dent Assoc. 2011
7. De Angelis F., Minnoni A., Vitalone L.M., Carluccio F., Vadini M., Paolantonio M., D'Arcangelo C.; *Bond strength evaluation of three self-adhesive luting systems used for cementing composite and porcelain*; University of Chieti, Department of Operative Dentistry, Dental School, Chieti, Italy Oper Dent. 2011 Nov-Dec
8. Peutzfeldt A., Sahafi A., Flury S.; *Bonding of restorative materials to dentin with various luting agents*; Department of Preventive, Restorative and Pediatric Dentistry, University of Bern, Bern, Switzerland. Oper Dent. 2011 May-Jun
9. Yüksel E., Zaimođlu A.; *Influence of marginal fit and cement types on microleakage of all-ceramic crown systems*; Department of Prosthodontics, Faculty of Dentistry, Ankara University, Ankara, Turkey; Braz. oral res. vol.25 no.3 São Paulo May/June 2011
10. Manso A.P., Silva N.R., Bonfante E.A., Pegoraro T.A., Dias R.A., Carvalho R.M.; *Cements and adhesives for all-ceramic restorations*; Department of Operative Dentistry, University of Florida, College of Dentistry; Dent Clin North Am. 2011 Apr

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Peruzzo Industrie Grafiche S.p.A.
Sweden & Martina S.p.A.